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CE 2777 (PPER(EU) 2016/425)

SATRA Technology Europe Ltd.
Bracetown Business Park
Clonee, D15 YN2P, Irlandia

INSTRUCTION FOR USE

POWDER-FREE NITRILE PROTECTIVE GLOVES

REFERENCE NUMBER

BECARE ART NITRILE COATS 504050**

BECARE ART NITRILE COATS GOODPACK+ 5046** / 504060**

**= 01 (size XS), 02 (size S), 03 (size M), 04 (size L), 05 (size XL)

	XS	S	M	L	XL
Available sizes	5/6	6/7	7/8	8/9	9/10

Read the Instruction for use before use

Version 01/2025

PURPOSE

Gloves for protection against chemicals and pathogens. Protective gloves for high risks. In addition, they have a positive effect on moisturising the wearer's skin.

INSTRUCTIONS FOR USE OF THE PRODUCT

Before use, after donning and during use, check the gloves for any defects or imperfections and discontinue use immediately if any damage, i.e. tears, abrasions, appears. Dry your hands before putting them on. Make sure that chemicals or chemical residues do not get through the sleeve. Hold the glove by the end with one hand when putting it on. Align the thumb of the glove with the thumb of your hand and then slide the other fingers into position. Pull the glove over the palm area to check the fit. Then do the same with the other glove. Removal: hold the glove by the end and pull towards the fingers until the glove is removed. Where hand hygiene is indicated prior to contact requiring gloves, wipe or wash hands before donning and removing gloves.

DISPOSAL

Used gloves may be contaminated and must be disposed of according to hospital and/or local regulations.

STORAGE

The shelf life of products stored as recommended is indicated on each package.

INCIDENT REPORTING

In the event of a serious incident related to the use of the product, the manufacturer and the competent authority responsible for supervision should be contacted. For further information, please contact the manufacturer, whose details are indicated at the end of this document.

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Measured breakthrough time in min.	>10	>30	>60	>120	>240	>480
Protection level	1	2	3	4	5	6

WARNINGS, PRECAUTIONS AND CONTRAINDICATIONS

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals and other factors affecting performance, such as temperature, abrasion, degradation, etc. Resistance to chemicals has been assessed under laboratory conditions on hand samples (except for gloves with a length equal to or longer than 400 mm - for

which the glove is also tested) and only applies to the chemicals tested. The result may vary if the chemicals are used in a mixture. It is advisable to check that the gloves are suitable for the intended use, as workplace conditions may differ from the tests carried out depending on temperature, abrasion and degradation. During use, protective gloves may provide less resistance to hazardous chemicals due to changes in physical properties. Movement, snagging, friction, degradation due to contact with the chemical, etc. can significantly reduce the actual service life. In the case of corrosive chemicals, degradation may be the most important factor to consider when selecting chemical-resistant gloves. The level of degradation (EN ISO 374-4:2019) indicates the change in puncture resistance of the glove after exposure to the chemical. The penetration resistance has been estimated under laboratory conditions and applies only to the test sample. Do not use if the glove is visibly torn or damaged. Gloves should be changed after each patient. Always choose the correct glove size to fit your hand size. For single use only. Poorly fitting gloves significantly reduce dexterity and cause hand fatigue. Using the wrong size glove leads to inadequate hand protection. Ingredients used in the manufacture of gloves may cause allergic reactions in some users. Consult a doctor immediately if allergic reactions occur. If required, a list of substances used in the manufacture of gloves that may cause allergic reactions is listed in Annex G of EN ISO 21420:2020, can be provided upon request.

EXPLANATION OF SYMBOLS AND PICTOGRAMS USED

Protective gloves against hazardous chemicals and micro-organisms Part 1: Terminology and chemical risk requirements. EN ISO 374-1. Definition of glove palm penetration time (1 $\mu\text{g}/\text{cm}^2/\text{min}$). Type A: level min. 2 for 6 chemicals, Type B level min. 2 for 3 chemicals, Type C level min. 1 for 1 chemical (no code under pictogram).

EXPLANATION OF THE MEANING OF THE PICTOGRAMS USED

 Permissible temperature	 Protect from light	 Store in a dry place	 Protect from ozone	 Suitable for food contact
EN 455-1  AQL level	 Contains potential type IV chemical allergens	 Do not reuse	autolayering  Autolayering: mechanical layering, alternating glove	 Does not contain natural rubber latex

EN ISO 374-5:2016  Virus Protective gloves against viruses, bacteria and fungi according to EN ISO 374-5:2016	EN ISO 374-1:2016+A1:2018/Type B  KPT Protective gloves against chemicals according to the standard EN ISO 374-1:2016 +A1:2018, Typ B	 Production date	 Manufacturer	 Importer
 Expiry date	 Goodpack layering technology	 User manual	 Dispose of in accordance with the regulations of the facility or the local regulations in force	 Indicator of where to open the packaging
 Code UDI	 Reference number	 Lot number	 Authorised representative in the EU	 Medical device in compliance with the Regulation UE 2017/745
 Packaging for recycling corrugated cardboard	 Packaging for recycling corrugated cardboard			

CONFORMITY ASSESSMENT

Latex-free: Yes; Powder-free: Yes

This product is a Category III PPE in accordance with Regulation (EU) 2016/425 and meets the standards: EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN 16523-1:2015+A1:2018, EN ISO 374-4:2019, EN ISO 374-5:2016, ISO 16604:2004.

Notified Body responsible for EU type examination (**Module B**) and conformity to type based on quality assurance of the production process (**Module D**): SATRA Technology Europe Ltd. (NB nr 2777), Bracetown Business Park, Clonee, D15YN2P, Ireland.

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SKAMEX S.A.
93-121, Łódź
Częstochowska 38/52
Poland

CE 2777 (PPER(EU) 2016/425)

SATRA Technology Europe Ltd.
Bracetown Business Park
Clonee, D15 YN2P, Ireland

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MDSS GmbH